



HOVERMATT® POSITIONING KIT

The HoverMatt® Positioning Kit includes the HoverMatt® SPU, Positioning Wedges and HoverCover™ Absorbent Disposable Cover, providing caregivers with all the tools needed to safely transfer, boost, position and offload patients without lifting or straining. The HoverMatt reduces the force required to move a patient by 80-90%, making patient handling easier and safer, while the Positioning Wedges effectively offload the sacrum and support the patient. The highly breathable HoverMatt SPU and HoverCover minimize heat and moisture buildup, making the patient more comfortable and protecting the integrity of the skin for improved all-day patient care.

HOVERMATT® POSITIONING KIT*



HoverMatt® SPU Benefits:

- Inflated mattress cradles the patient while moving them in a stable position, eliminating skin shear and bruising and providing maximum comfort to the patient
- Air-assisted technology is proven to reduce workplace injuries and associated costs by 50-100% when implemented into a Safe Patient Handling Program

Product Specifications

MATERIAL SPECS: Non-woven polypropylene fiber, Nylon fabric, Latex free

DIMENSIONS: 34" W x 78" L; 39" W x 78" L; 50" W x 78" L;

34" W x 45" L; 39" W x 45" L


Positioning Wedge Benefits:

- High resiliency foam supports the patient, while minimizing additional pressure from compression over time.
- 30° position provides effective sacral offloading
- Non-slip cover secures and stabilizes the device under the patient
- Patient dedicated device reduces potential for cross-contamination
- Polyurethane and nylon cover simplifies cleaning

Product Specifications

MATERIAL SPECS: Polyurethane foam, Latex free

DIMENSIONS: 15" L x 14" W x 7" H, 30°



HoverCover™ Benefits:

- Breathable cover wicks away moisture to help control the microclimate of the patient's skin
- Highly absorbent layer protects the HoverMatt from fluids to help with infection control and facilitate cleanup

Product Specifications

MATERIAL SPECS: Non-Woven, Breathable PE, Latex free

DIMENSIONS: 65" L x 39" W



Kit Model Numbers:

HPK-34SPU

HPK-39SPU

HPK-50SPU

HPK-34SPU-HLF

HPK-39SPU-HLF



4482 Innovation Way
Allentown, PA 18109

800.471.2776

Fax 610.694.9601

www.HoverMatt.com

Info@HoverMatt.com

* Each kit contains (1) HoverMatt SPU, (2) Positioning Wedges,
and one pack of (5) HoverCovers.

BRO_HPK_1117