

HOVERCOVER™

Breathable • Ultra Absorbent • Protective



ELEVATE YOUR STANDARD OF CARE

Elevating Staff Safety, Patient Care, and Value

Visit HoverMatt.com for More Information




HOVERCOVER™

ULTRA ABSORBENT • INCONTINENCE CARE • SUPERIOR PROTECTION

HoverCover Benefits:

- Breathable with high absorbency
- Draws away moisture to keep patient cool and dry
- Stays in place to reduce leaks and wrinkles
- Compatible with low air loss mattresses
- Larger size prevents overlapping of multiple chux
- Single-patient use, full and half size



Superior protection for air-assisted products:
 HoverMatt® SPU, HoverSling®, HoverMatt® Half-Matt,
 & Q2Roller® Lateral Turning Device.

Model #: HMAc-100-5PACK (5 ea.)
 Model #: HMAc-100-CASE (50 ea.)
 39" W x 65" L

Model #: HTR-CHX-5PACK (5 ea.)
 Model #: HTR-CHX-CASE (50 ea.)
 40" W x 40" L

One full-length HoverCover vs. three average-sized chux

	HoverCover	Average chux*
1 Size (per chux)	39" x 65"	30" x 36"
2 Quantity needed to cover bed	1	3
3 Absorbency (per chux)	6.7 L	1.5 L
4 Reduced wrinkles/movement	✓	✗
5 No leaking between overlap	✓	✗
6 One product for quick bed set up	✓	✗
7 Keeps surface clean under product (bed, air transfer mattress, sling, etc.)	✓	✗

*Ultra Absorbency