



HOVERSLING[®]

Repositioning Sheet

The breathable HoverSling[®] Repositioning Sheet functions as both a HoverMatt SPU and sling to streamline patient handling tasks and reduce the amount of time and money spent using multiple products. With the HoverSling Repositioning Sheet, caregivers can use one product for lateral transfers, boosting, positioning, turning and vertical lifting. The patient also benefits from trusted HoverMatt air technology that reduces skin shear during transfers and repositioning.

HOVERSLING®

Repositioning Sheet



HoverSling® Repositioning Sheet Benefits

- Highly breathable fabric allows for greater moisture vapor exchange and quicker moisture evaporation
- Independent laboratory testing proves compatibility with low air loss surfaces. Can be left under patient for all day care
- Combination HoverMatt and Repositioning Sheet reduces inventory and offers streamlined approach for patient transfers and positioning
- Supine sling design fits most patients, simplifying the decision-making process for caregivers



TURN



TRANSFER



LIFT

Product Specifications

MATERIAL SPECS

Sewn construction

Top: Breathable polyester knit fabric

Bottom: Breathable nylon fabric

Latex free

DIMENSIONS

39" W x 78" L – HMSLING-39RS-B (5 per box)

50" W x 78" L – HMSLING-50RS-B (5 per box)

CAN BE USED WITH ALL LOOP STYLE HANGER BARS

Weight limit 1,000 pounds

Patent No(s). 6,898,809; 7,712,170; 7,373,680



4482 Innovation Way
Allentown, PA 18109

800.471.2776
Fax 610.694.9601

www.HoverMatt.com
Info@HoverMatt.com

BRO_RS_0621