



HOVERMATT®
SPU

HOVERSLING®

Lateral Transfers • Boosting • Turning • Positioning • Proning • Lifting

FACILITATE ALL-DAY PATIENT CARE WITH THE BREATHABLE HOVERMATT® SINGLE-PATIENT USE AIR TRANSFER SYSTEM.

Designed to be multifunctional, the breathable HoverMatt® Single-Patient Use (SPU) Air Transfer System makes patient transfers, boosting, positioning and proning easier, while protecting caregiver safety. A cushion of air beneath the inflated HoverMatt SPU reduces the force required to move a patient by 80-90%, enabling caregivers to safely transfer patients without lifting or straining. By reducing injuries related to lateral transfers and repositioning, the HoverMatt SPU helps to improve staff retention while meeting legislative guidelines for safe patient handling.

Improves Patient Comfort & Safety

- Advanced breathable textiles enable caregivers to confidently leave the product under their patients
- Inflated mattress cradles the patient while moving them in a stable position, eliminating skin shear and bruising and providing comfort to patients in pain

Reduces Cost

- Single-patient use product eliminates costs and logistics associated with laundering, while greatly reducing the risk of cross-contamination
- Air-assisted technology is proven to reduce workplace injuries and associated costs 50-100% when implemented into a Safe Patient Handling Program

Facilitates Care Throughout the Facility

- With a 1200 lb. weight limit, one device can be used for multiple patient handling tasks across the facility, including in Radiology, Diagnostic Imaging and ancillary departments
- Supports Safe Patient Handling recommendations and guidelines for increased healthcare worker safety from OSHA, AORN, NAON and ANA

HT-Air® 1200 Air Supply offers six airflow options, including a setting that can be used to inflate the HoverMatt slowly for sensitive patients and ensure proper positioning prior to full inflation.



HOVERMATT®
SPU

HoverMatt SPU with unique connecting straps

The breathable HoverMatt® SPU Link features all of the trusted benefits of the original HoverMatt SPU combined with unique connecting straps that attach to the bed helping to keep the device in place to reduce the need for boosting and improve patient comfort.



HOVERMATT®
SPU Link

Ortho, L&D and UroGyn

Specifically designed for use in the OR and Labor and Delivery, the HoverMatt® SPU Half-Matt delivers the same safe and easy patient transfers as the standard HoverMatt but in a shorter length. At 45" long, the SPU Half-Matt ends right below the buttocks to accommodate specialty tables. It is ideal for use during arthroscopy, cystoscopy, labor and delivery, and other gynecological or urological procedures as an effective yet unobtrusive patient transfer device.



HOVERMATT®
SPU Half-Matt

Bariatrics and UroGyn

The HoverMatt® SPU Split-Leg Matt is a specialty transfer mattress designed for use in the lithotomy position or on surgical tables when split-leg attachments are used. The SPU Split-Leg Matt has a patented design that separates the bottom half of the transfer mattress into two individual segments, or "legs." These legs snap together for effortless lateral patient transfers and are then unsnapped for convenient positioning on the procedure table.



HOVERMATT®
SPU Split-Leg

HOVERTECH ADVANTAGES

- Innovative market leader focused on customer needs
- Extensive air-assisted breathable transfer device product line
- Comprehensive compatibility safety testing:
 - Radiolucency
 - Pressure Injury Prevention
 - Electrosurgery Compatibility
 - Static Testing
 - MRI Compatibility
 - Attenuation
 - OR Particle Test
 - Heat Transfer
 - Flammability
- Equipment acquisition programs including the HT-Air Placement Program, HT-Air Rental Program, and HoverTech Rewards Program
- Professionally trained, account-dedicated sales and customer service representatives
- Safe Patient Handling Healthcare Ergonomics Lifting Program (HELPSM) provides consultative services to assist facilities to develop or improve a SPHM program

STREAMLINE PATIENT HANDLING TASKS

The breathable HoverSling® Repositioning Sheet functions as both a HoverMatt SPU and sling to streamline patient handling tasks and reduce the amount of time and money spent using multiple products. With the HoverSling Repositioning Sheet, caregivers can use one product for lateral transfers, boosting, positioning, turning and vertical lifting. The patient also benefits from trusted HoverMatt air technology that reduces skin shear during transfers and repositioning.



HOVERSLING®
Repositioning Sheet

The breathable HoverSling Split-Leg offers all of the benefits of the HoverMatt SPU Split-Leg for lateral transfers, boosting, turning and repositioning. When used as a lift sling for seated transfers, the patented design cradles the patient's limbs evenly and securely for greater patient comfort. This unique product can also be used as a limb lifter or stabilizer for wound care and Labor and Delivery applications.



HOVERSLING®
Split-Leg

HOVERMATT® PRODUCT SPECIFICATIONS

Sewn construction

Top: Non-woven polypropylene fiber

Bottom: Nylon fabric

Latex free

HoverMatt Single-Patient Use

34" W x 78" L – HM34SPU-B (10 Per Box)

39" W x 78" L – HM39SPU-B (10 Per Box)

50" W x 78" L – HM50SPU-B (5 Per box)

HM50SPU-1Matt-B (1 Unit)

Weight limit 1,200 pounds

Patent No(s). 7,373,680; 7,712,170

HoverMatt SPU Link

34" W x 78" L – HM34SPU-LNK-B (10 Per Box)

39" W x 78" L – HM39SPU-LNK-B (10 Per Box)

Weight limit 1,200 pounds

Patent No(s). 6,898,809; 7,712,170; 7,373,680

and patents pending

HoverMatt SPU Half-Matt

34" W x 45" L – HM34SPU-HLF-B (10 Per Box)

39" W x 45" L – HM39SPU-HLF-B (10 Per Box)

Weight limit 600 pounds

Patent No. 6,073,291

HoverMatt SPU Split-Leg Matt

34" W x 70" L – HM34SPU-SPLIT-B (10 Per Box)

Weight limit 1,200 pounds

Patent No(s). 7,266,852; 7,406,723; 7,574,761

HOVERSLING® PRODUCT SPECIFICATIONS

Sewn construction

Top: Polyester knit

Bottom: Nylon fabric

Latex free

HoverSling Split-Leg

34" W x 70" L – HMSLING-34-B (10 Per box)

Weight limit 700 pounds

Patent No(s). 7,373,680; 8,566,977; 7,266,852;

7,574,761; 7,712,170

HoverSling Repositioning Sheet

39" W x 78" L – HMSLING-39RS-B (5 Per Box)

50" W x 78" L – HMSLING-50RS-B (5 Per Box)

Weight limit 1,000 pounds

Patent No(s). 6,898,809; 7,712,170; 7,373,680

and patents pending

Compatible with HoverTech Air Supplies

Air200G (800W), Air400G (1100W), and HTAIR1200/2300

(Variable Speed)

If the devices are laundered, the "DO NOT WASH" label at the foot end of the product will dissolve and state "DISCARD IMMEDIATELY".



4482 Innovation Way
Allentown, PA 18109

800.471.2776
Fax 610.694.9601

www.HoverMatt.com
Info@HoverMatt.com

BRO_SPU_0621