

SAFE PATIENT HANDLING

Minimum Crew.

- No waiting for additional personnel or equipment when every second counts
- Increase completed calls per shift and generate more revenue due to faster transport times



Maximum Safety.

- Lift, transfer, and transport patients without risking the safety of your crew or patient
- Minimize manual lifting and decrease costly workers' compensation injuries



Evacuation EMS HoverJack® Device

The Evacuation EMS HoverJack® Device is HoverTech's all-in-one transport, lift, and stair evacuation solution that minimizes many of the challenges EMS crews face in taking patients from their home to the hospital.

Evacuation HoverJack® Device Benefits:

- Air chambers inflate to lift the patient from the floor in a supine or seated position, creating a no-lift environment for the crew
- Tapered head-end, adjustable inflation chambers, and steering straps improve maneuverability around corners and down most stairwells
- Quick-release valves expedite deflation to lower patient when needed
- Support skirt with replaceable Teflon coated bottom allows device to be used across rough terrain
- Weight capacity of 700 lbs. for stairwell evacuations or 1200 lbs. for vertical lifts accommodates most patients
- Dimensions: 32" W x 78" L – HJ32EV-2



HoverJack® Battery Cart

Houses a rechargeable battery and onboard outlet to power an Air Supply.

HoverMatt® Air Transfer System

The HoverMatt® Air Transfer System allows EMTs to laterally transfer patients from their bed to the Evacuation EMS HoverJack or to a stretcher with ease. When inflated, a cushion of air beneath the mattress reduces the force needed to move a patient by 80–90%. This reduction in force minimizes crew and patient injuries that are often caused by lateral transfers.

HoverMatt® Air Transfer System Benefits:

- Air transfer technology allows first responders to laterally transfer patients quickly and safely, while reducing patient skin shear and bruising
- Weight limit of 1200 lbs. accommodates most patients
- Available in the following options and sizes:
 - Single-Patient Use: 34" W x 78" L – HM34SPU-B (10 per box); 39" W x 78" L – HM39SPU-B (10 per box); 50" W x 78" L – HM50SPU-B (5 per box)
 - Reusable: 34" W x 78" L – HM34HS (Each); 39" W x 78" L – HM39HS (Each); 50" W x 78" L – HM50HS (Each); 34" W x 78" L – HM34DC (Each)
 - Double-Coated: 39" W x 78" L – HM39DC (Each); 50" W x 78" L – HM50DC (Each)

