



HOVERJACK[®]

Lifting • Evacuation • Vehicle Extraction

FACILITATE PATIENT LIFTING AND EVACUATIONS IN THE MOST CHALLENGING ENVIRONMENTS WITH THE HOVERJACK® AIR PATIENT LIFT.

The HoverJack® Air Patient Lift allows caregivers to safely lift patients who have fallen without gathering a lift team. The HoverJack inflates to lift patients from the floor to bed or stretcher height in a supine position, maximizing patient comfort and minimizing the risk of injury to both the patient and caregivers. Weight capacity: 1200 lbs. The HoverJack is a critical component for fall prevention and emergency preparedness.

Promotes Patient Comfort & Dignity

- Supine lifting feature prevents patient from incurring further injuries after a fall
- Air inflated chambers provide a more comfortable lift compared to other methods
- Available in two widths with a weight capacity of 1200 lbs. to accommodate most patients

Improves Staff Efficiency & Safety

- Air-assisted technology reduces the number of caregivers required to lift a patient when responding to a fall
- Self-sealing inflation valves facilitate ease of use, while the deflation valves allow users to quickly lower the patient if needed

Facilitates Care In & Out of the Facility

- System flexibility allows for use anywhere in the facility including small spaces, like the bathroom, where floor lifts or ceiling lifts cannot be utilized
- Individual chamber inflation design allows users to easily transfer patients to surfaces of varying heights
- Teflon® coated bottom facilitates movement across all surfaces, including parking lots

The HoverMatt® Air Transfer System can be used with the same air supply for safe and easy lateral patient transfers from the HoverJack to a bed or stretcher.



HOVERJACK®

Simplify Stairwell Evacuations with the Evacuation HoverJack® Device

The HoverJack® Evacuation Lift allows caregivers to perform safe and efficient stairwell evacuations while keeping the patient in a supine position. Adjustable air chambers inflate to provide a stable and comfortable surface for patients, while easily navigating doorways and narrow spaces. The Teflon® coated bottom allows caregivers to effectively transport the HoverJack Evacuation device across most surfaces and down stairwells. With all of the same features as the standard HoverJack, the HoverJack Evacuation has the following added benefits to accommodate emergency evacuations.

- Cross-over chest straps, dual lateral body straps and foot-end pouch provide patient security during evacuation
- Tapered head-end improves maneuverability around corners and down all stairwell designs, including curving, L-shaped and steep, straight stairs
- Inflatable head-end wedge facilitates adjustable positioning to improve patient's breathing and comfort
- Weight capacity of 700 lbs. for stairwell evacuations and 1200 lbs. for vertical lifts
- Personal belongings compartment built into device offers convenient storage



HOVERJACK®
Evacuation



The HoverJack may also be used in combination with the HoverMatt to safely and easily extract a patient from a vehicle.



HT-Air® 1200 Air Supply

The HT-Air® 1200 Air Supply offers two air flow speeds and a variable speed air flow option. The HT-Air can be used with our full range of products to maximize product efficiency.

HoverJack® Cart

The HoverJack® Cart is a rolling storage system that allows caregivers to easily respond in the event of a patient fall. With a standard size top basket and larger bottom basket, this cart holds the HoverJack, Air Supply and a HoverMatt transfer mattress.



HoverJack® Battery Cart

The HoverJack® Battery Transport System offers a durable, ergonomic solution for portability and storage of the Air200G/Air400G Air Supply and any HoverJack device, featuring the added convenience of a built-in, rechargeable battery.

With its spacious, lightweight and easy to steer design, the HoverJack Battery Transport System is an all-in-one resource for Rapid Response Teams.

PRODUCT SPECIFICATIONS

HoverJack Air Patient Lift

Topside Material: Nylon Oxford
Underside Material: Teflon® infused polyester
Latex free

AVAILABLE SIZES

32" W x 72" L – HJ3202T

39" W x 77" L – HJ3902T

Compatible with all HoverTech Air Supplies

Evacuation HoverJack

Topside Material: Nylon Oxford
Foot Pocket: Urethane coated Packcloth
Bottom Chamber: Reinforced with urethane coated Packcloth
Underside Material: Teflon® infused polyester

Size: 32" W x 72" L – HJ32EV

Compatible with all HoverTech Air Supplies

Carts

HoverJack Cart – HJC-200

HoverJack Battery Cart – HJBSC-200*

**Please note this item is not compatible with the HTAIR1200 air supply*



4482 Innovation Way • Allentown, PA 18109

Phone: 800.471.2776 • Fax: 610.694.9601 • info@HoverMatt.com
HoverMatt.com