



# PRODUCT CATALOG

# ELEVATE YOUR STANDARD OF CARE

Elevating staff safety, patient care, and value

## HOVERMATT®

The HoverMatt® Reusable Air Transfer Mattress makes patient transfers easier for caregivers and patients. When inflated, a cushion of air beneath the mat reduces the pull force required to move the patient by 80 – 90%, enabling caregivers to transfer, turn and position patients safely. A Double-Coated model is ideal for use when excessive bodily fluids are a concern. Weight capacity: 1200 lbs.

### Reusable Model Numbers

34" W x 78" L – HM34HS (Each)

39" W x 78" L – HM39HS (Each)

50" W x 78" L – HM50HS (Each)

34" W x 78" L – HM34DC (Each) *Double-Coated*

39" W x 78" L – HM39DC (Each) *Double-Coated*

50" W x 78" L – HM50DC (Each) *Double-Coated*



## HOVERMATT® SPU

The HoverMatt® Single-Patient Use Air Transfer Mattress works just like the HoverMatt Reusable, using a cushion of air that enables caregivers to safely and easily transfer, turn and position patients. Single-Patient Use (SPU) products are made with highly breathable fabric, making them safe and beneficial to leave under patients for all-day care that supports pressure injury prevention. Weight capacity: 1200 lbs.

### Model Numbers

34" W x 78" L – HM34SPU-B (10 per box)

39" W x 78" L – HM39SPU-B (10 per box)

50" W x 78" L – HM50SPU-B (5 per box)



## HOVERMATT® SPU Link

The HoverMatt® Single-Patient Use Link Air Transfer Mattress features all of the trusted benefits of the original HoverMatt SPU combined with unique connecting straps that attach to the bedframe. The connecting straps reduce the need for boosting and improve patient comfort, and serve as a tool to facilitate patient egress. Weight capacity: 1200 lbs.

### Model Numbers

34" W x 78" L – HM34SPU-LNK-B (10 per box)

39" W x 78" L – HM39SPU-LNK-B (10 per box)

50" W x 78" L – HM50SPU-LNK-B (5 per box)





# HOVERMATT®

## Half-Matt

The HoverMatt® Single-Patient Use Half-Matt and Reusable Half-Matt are both 45" in length to provide an effective yet unobtrusive solution for specialty surgical tables. These mats are ideal for use during cystoscopy, labor and delivery and other gynecological, urological and orthopedic procedures. Weight capacity: 600 lbs.

### Single-Patient Use Model Numbers

34" W x 45" L – HM34SPU-HLF-B (10 per box)

39" W x 45" L – HM39SPU-HLF-B (10 per box)

### Reusable Model Numbers

34" W x 45" L – HM34HS-MINI (Each)

34" W x 45" L – HM34DC-MINI (Each) *Double-Coated*



# HOVERMATT®

## SPU Split-Leg

The HoverMatt® SPU Split-Leg Matt is a specialty transfer, turn, and positioning mat designed for use in the lithotomy position or on surgical tables when split-leg plates are employed. Weight capacity: 1200 lbs.

### Model Number

34" W x 70" L – HM34SPU-SPLIT-B (10 per box)



# HOVERSLING®

The breathable HoverSlings function as both a HoverMatt SPU and lift sling to streamline patient handling tasks and reduce the time and money spent using multiple products. In addition to safely enabling lateral transfers and turning and positioning patients, the HoverSling® Repositioning Sheet is also designed for vertical supine lifting, while the HoverSling® Split-Leg Device is designed for vertical seated lifting.

Weight capacities: 700 - 1000 lbs. See below for details.

### HoverSling Repositioning Sheet Model Numbers

39" W x 78" L – HMSLING-39RS-B (5 per box) - Weight Capacity: 1000 lbs.

50" W x 78" L – HMSLING-50RS-B (5 per box) - Weight Capacity: 1000 lbs.

### HoverSling Split-Leg Model Numbers

HMSLING-34-B – 34" W x 70" L (10 per box) - Weight Capacity: 700 lbs.



# HOVERJACK®

The HoverJack® Air Patient Lift allows caregivers to safely lift patients who have fallen without gathering a lift team. The HoverJack inflates to lift patients from the floor to bed or stretcher height in a supine position, maximizing patient comfort and minimizing the risk of injury to both the patient and caregivers. Weight capacity: 1200 lbs.

## Model Numbers

32" W x 72" L – HJ3202T (Each)  
39" W x 77" L – HJ3902T (Each)



# HOVERJACK®

## Evacuation

The HoverJack® Evacuation Lift is an emergency stairwell evacuation solution designed for patients who cannot be evacuated using a bed mattress. The specially designed safety straps and foot-end pouch hold the patient securely, while the Teflon-coated bottom allows for easy transport across most surfaces and down stairs. Weight capacity: 700 lbs. for stairwell evacuations and 1200 lbs. for vertical lifts.

## Model Number

32" W x 72" L – HJ32EV (Each)



# HOVERJACK®

## EMS Evacuation

The HoverJack® EMS Evacuation Lift is an air-assisted transport, lift, and stair evacuation solution that minimizes the challenges EMS crews face in transporting bariatric patients. The reinforced, replaceable bottom chamber skirt and harness-style safety straps make this version appropriate for the heavy demands of the EMS market. Weight capacity: 700 lbs. for stairwell evacuations and 1200 lbs. for vertical lifts.

## Model Number

32" W x 72" L – HJ32EV-2 (Each)



# Q2ROLLER®

The Q2Roller® Lateral Turning Device comprises two breathable inflatable chambers that facilitate easy and efficient turning for patient care and pressure relief on bony prominences. By eliminating the need for manual lifting, the Q2Roller protects both the staff from injury and the patient from friction and shearing injuries associated with traditional positioning methods. Weight capacity: 600 lbs.

## Dimensions

Deflated Chambers: 44.5" L x 46.5" W (total)  
Apron: 44.5" L x 72" W

## Model Number

HTR-200 (5 per box)  
Each Q2Roller comes with one pack of 5 Q2R Pads:  
(HTR-CHX-5PACK)





# ACCESSORIES

## AIR SUPPLIES

### HT-AIR® 1200

The HT-Air® 1200 Air Supply offers two air flow speeds and a variable speed air flow option. The HT-Air can be used with our full range of products to maximize product efficiency.

**Model Number** HTAIR1200 (Each) *Not compatible with HoverJack Battery Cart*



### AIR400G

The Air400G Air Supply offers rapid inflation of HoverMatt and HoverJack products. Compatible with the HoverJack Battery Cart.

**Model Number** AIR400G (Each)



### MRI HOSE

The 25' MRI Hose enables caregivers to use the HoverMatt in MRI by keeping the air supply out of the imaging environment.

**Model Numbers** HTA-MRI (Each) – for HTAIR1200 AS-MRI (Each) – for Air400G



### AIR SUPPLY HOSE COVERS

The Air Supply Hose Cover is a convenient, disposable cover that protects the entire length of the Air Supply Hose from soiling. It is especially useful for infection control and high fluid areas.

**Model Number** AS-HC (25 per package)



## POSITIONING WEDGES

The 30 degree Positioning Wedges comply with NPIAP Guidelines for the prevention and treatment of pressure injuries. They feature antimicrobial, hypoallergenic and latex-free high density foam.

### Dimensions

15"L x 11"W x 7"H, 30°  
15"L x 11"W x 7"H, 30°  
20"L x 11"W x 7"H, 30°  
26"L x 12"W x 8"H, 30°

### Model Numbers

FPW-100 (2 per box) *Single-patient use*  
FPW-R-15S (Each) *Reusable*  
FPW-R-20S (Each) *Reusable*  
FPW-RB-26S (Each) *Reusable*

### Weight Capacity

Up to 350 lbs.  
Up to 350 lbs.  
Up to 350 lbs.  
Bariatric - 300 to 650 lbs.



# CARTS

## HOVERMATT® UNIVERSAL CART

The HoverMatt® Universal Cart is designed to hold all HoverTech air supply models and one box of HoverMatt SPU's or several HoverMatt Reusables.

Model Number HMC-U (Each)



## HOVERJACK® CART

The HoverJack® Cart is designed to hold the HoverJack Air Patient Lift, a HoverTech air supply and a HoverMatt to allow caregivers to respond with the right equipment in the event of a patient, staff or visitor fall.

Model Number HJC-200 (Each)



## HOVERJACK® BATTERY CART

The HoverJack® Battery Cart houses a rechargeable battery and onboard outlet to power an Air Supply (for use with Air200G and Air400G HoverTech Air Supplies only).

Model Number HJBSC-200 (Each)



## HOVERCOVER™ BREATHABLE ABSORBENT PAD

The HoverCover™ Breathable Absorbent Pad provides an absorbent and breathable layer that can be placed on top of any HoverMatt, HoverSling, or other patient handling products to protect them from becoming soiled. Available in Full and Half sizes.

### Model Numbers

#### FULL

HMAC-100-5PACK (5 per package)

HMAC-100-CASE (50 per case)

#### HALF

HTR-CHX-5PACK (5 per package)

HTR-CHX-CASE (50 per case)



## DISPOSABLE SHEET

The Disposable Sheet is electronically weakened down the center to allow for easy removal. It can be placed on top of the HoverMatt, as well as other patient handling products, to protect them from becoming soiled. It features an absorbent top and impervious underside material.

### Model Numbers

HMDS3-5PACK (5 per package)

HMDS3-CASE (100 per case)





## Maximize the Return on Investment of Your Safe Patient Handling Equipment

HoverTech offers the HELP<sup>SM</sup> Healthcare Ergonomic Lifting Program, which provides consultative services to assist with developing and/or improving a Safe Patient Handling and Mobility (SPHM) program.

- Clinical consulting experts work with your facility to identify and reduce risk factors associated with patient handling tasks, thereby minimizing workplace injuries and lowering workers' compensation and associated costs.
- The systematic implementation process promotes hospital ownership and skill development within your facility.
- The program is both flexible and scalable to meet your needs, whether starting a new SPHM program or enhancing an existing one.



4482 Innovation Way • Allentown, PA 18109

Phone: 800.471.2776 • Fax: 610.694.9601 • [info@HoverMatt.com](mailto:info@HoverMatt.com)

[HoverMatt.com](http://HoverMatt.com)

BRO\_CATALOG\_0821