



# HOVERMATT® T-BURG™

Trendelenburg Stabilization & Air Transfer Mattress

Patient Safety. Your Safety. One Device.



## ELEVATE YOUR STANDARD OF CARE

Elevating Staff Safety, Patient Care, and Value

# ONE PRODUCT FROM THE O.R. THROUGH RECOVERY

## Patient Stabilization. Patient Comfort. Staff Safety.

The HoverMatt® T-Burg™ is designed for safe Trendelenburg positioning *and* operates as an air-assisted device for safe patient transfers. The HoverMatt T-Burg supports the entire patient journey before, during, and after surgery.

### During Surgery

- **Breathable HoldFast™ fabric prevents** unwanted patient movement in steep Trendelenburg.
- **Efficiently** attaches to the O.R. Table in **< 60 seconds**.
- **HoldFast™** supports a microclimate that contributes to skin integrity while remaining under the patient for their entire stay.
- Safe & Efficient **patient air-transfer technology** from the O.R. to the nursing floor.



**What makes a great Trendelenburg Stabilization Mat?**

**No Unwanted Patient Movement**

**Less than 60 seconds**  
to efficiently set up

**Increased Staff Safety**  
when transferring a patient from the O.R. table to a bed or stretcher

### After Surgery

- **Gentler** on the patient.
- **No unnecessary** patient movement compared to other stabilization pads.
- **No patient “log rolling”** required after surgery to move from the O.R. table to recovery.
- **Placement of a pad or chux** after surgery is optional.



# PATIENT POSITIONING IN STEEP TRENDELENBURG & MORE

## Trendelenburg Positioning

With or without Lithotomy




## Reverse Trendelenburg Positioning



## Lateral Transfer

From OR table to a bed or stretcher



Patient Safety. Your Safety. One Device.

# THE LEADERS IN SAFE PATIENT HANDLING

Extensive offering of air-assisted device solutions for a broad range of patient handling and pressure injury prevention needs.

## HOVERMATT® AIR TRANSFER SYSTEM



HOVERMATT™  
T-Burg™



HOVERMATT™  
SPU



HOVERMATT™  
SPU Link



HOVERMATT™  
Half-Matt



HOVERMATT™  
SPU Split-Leg



HOVERMATT™

## HOVERSLING®



HOVERSLING™  
Repositioning Sheet



HOVERSLING™  
Split-Leg

## HOVERJACK®



HOVERJACK™



HOVERJACK™  
Evacuation

## Q2ROLLER® LATERAL TURNING DEVICE



Q2ROLLER™



WANT TO LEARN MORE ABOUT  
THE HOVERMATT® T-BURG™?

Scan the QR Code to visit our website

HOVERMATT™  
T-Burg™

Model Numbers:

HM-TBURG-34: 5 PC / PE BOX

HM-TBURG-34-KIT: Kit includes 5 ea

T-Burg Matts + 5 pair arm protectors

Weight Limit - 400 lbs. from 0 to 40°

Compatible with All HoverTech Air Supplies



4482 Innovation Way • Allentown, PA 18109

Phone: (800) 471-2776 • Fax: 610-694-9601 • [info@hovermatt.com](mailto:info@hovermatt.com)

[HoverMatt.com](http://HoverMatt.com)

Bro\_TBurg\_030123